Control FLQ Vitamin D**REF** K082M-10

10 x 2 mL

LOT 07045Aalto Scientific Ltd
230 Technology Pkwy
Eatonton, GA 31024
USAMDSS GmbH
Schiffgraben 41
30175 Hanover
Germany

2027-02-26

ENGLISH**INTENDED USE**

Control FLQ Vitamin D is an assayed, stable, bi-level reference control for use with assays designed to quantitate Vitamin D. It is intended to simulate human patient samples for the purpose of monitoring the precision of laboratory testing procedures for Vitamin D assays. When used for quality control purposes, it is recommended that each laboratory establish its own means and acceptable ranges and use the values provided only as guides.

This product is intended for use with quantitative assays on the indicated analyzer provided in the labeling. The Control FLQ Vitamin D should not be used for calibration or standardization of the Vitamin D assay. The Control FLQ Vitamin D is for In Vitro Diagnostic use only.

SUMMARY AND PRINCIPLE

Good laboratory practices require that stable reference materials be used to verify the accuracy and precision of testing methods and techniques. Control FLQ Vitamin D may be used as one would use human serum to obtain the stated Vitamin D values. This control will assist in the evaluation of proper performance of Vitamin D assays.

REAGENTS

Control FLQ Vitamin D is prepared using human serum with purified Vitamin D added to achieve the desired concentration levels. Control FLQ Vitamin D is manufactured according to standard quality control procedures. The manufacturer guarantees stability and consistency of this product.

WARNINGS AND PRECAUTIONS

Because this product is of human origin, it should be handled as though capable of transmitting infectious diseases. Each serum, plasma or whole blood donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled at the Centers for Disease Control's Biosafety Level 2.

Control FLQ Vitamin D is intended solely for in vitro diagnostic use for the purpose described on the labeling. Audit® MicroControls™, Inc. shall not be liable for any unclaimed damages arising from any other usage.

MATERIALS PROVIDED

Control FLQ Vitamin D, 10 x 2 mL

STORAGE AND STABILITY

Control FLQ Vitamin D may arrive thawed. It should be stored at -15°C or colder upon receipt until use.

Control FLQ Vitamin D is stored at -15°C or colder and will remain stable in the unopened bottle until the expiration date. Do not store in a frost-free freezer. After thawing, the contents should be used according to the instrument manufacturer's instructions and stored at 2-8°C for 30 days.

PROCEDURE

Follow the manufacturer's instructions provided for Vitamin D procedures. Verify that the lot number on the vial matches the assay sheet. To avoid evaporation, do not leave the vial uncapped. Controls should be run:

1. daily, in conjunction with patient samples.
2. as recommended by the instrument manufacturer.
3. as required by the relevant regulatory agency.

INSTRUCTIONS FOR USE

1. Remove the bottles from the freezer and allow to thaw at room temperature. Do not shake. Do not mix mechanically.
2. Mix gently by inversion to assure complete mixing of the contents.
3. Remove cap and invert bottle. Gently squeeze the dropper bottle.
4. After sampling, wipe off dropper tip, replace cap, and store at 2-8°C.

EXPECTED VALUES

The performance range for each level, based on data by combining estimates of assay variance as determined by participating laboratories using approved FDA instruments and reagents, is provided below. Average values obtained in the laboratory should fall within the performance range although the recovery may not be identical with the mean value listed. Variation between labs will be greater than the precision for any one instrument. Accuracy and precision depend on differences in equipment, reagents, supplies and techniques. Therefore, a lab must establish its own acceptable target values and ranges.

LIMITATIONS OF THE PROCEDURE

The Control FLQ Vitamin D should not be used for calibration or standardization of the Vitamin D assays.

ORDERING INFORMATION

PRODUCT NUMBER	PRODUCT DESCRIPTION	PRODUCT PACKAGING
K082M-10	Control FLQ Vitamin D	10 x 2 mL

Distributed by AUDIT MicroControls™ Inc. - U.S. customers only please call (866) 252-8348 or www.auditmicro.com

Catalog Number

For In Vitro
Diagnostic UseUse By
(YYYY-MM-DD)

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