Control FD BNP

REF K023M-6

6 x 1 mL

LOT 07173

 2027-07-17

Aalto Scientific Ltd
230 Technology Pkwy
Eatonton, GA 31024
USA

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ENGLISH

INTENDED USE

Control FD BNP is an assayed, stable, lyophilized, bi-level reference control for use with assays designed to quantitate BNP and NT-proBNP. It is intended to simulate human patient samples for the purpose of monitoring the precision of laboratory testing procedures for BNP and NT-proBNP assays. When used for quality control purposes, it is recommended that each laboratory establish its own means and acceptable ranges and use the values provided only as guides.

This product is intended for use with quantitative assays on the indicated analyzers provided in the labeling. The Control FD BNP should not be used for calibration or standardization of the BNP and NT-proBNP assay. The Control FD BNP is for In Vitro Diagnostic use only.

SUMMARY AND PRINCIPLE

Good laboratory practices require that stable reference materials be used to verify the accuracy and precision of testing methods and techniques. Control FD BNP may be used as one would use human plasma to obtain the stated BNP and NT-proBNP values. This control will assist in the evaluation of proper performance of BNP and NT-proBNP assays.

REAGENTS

Control FD BNP is prepared using human plasma with purified BNP and NT-proBNP added to achieve the desired concentration levels, and 0.1% sodium azide as a preservative. Control FD BNP is a lyophilized product manufactured according to standard quality control procedures. The manufacturer guarantees stability and consistency of this product.

WARNINGS AND PRECAUTIONS

Because this product is of human origin, it should be handled as though capable of transmitting infectious diseases. Each serum, plasma or whole blood donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled at the Centers for Disease Control's Biosafety Level 2.

This product contains less than 0.1% sodium azide that may react with lead and copper plumbing to form potentially explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Control FD BNP is intended solely for in vitro diagnostic use for the purpose described on the labeling. Audit® MicroControls™, Inc. shall not be liable for any unclaimed damages arising from any other usage.

MATERIALS PROVIDED

Control FD BNP, 6 x 1 mL

STORAGE AND STABILITY

Control FD BNP is stored at 2-8°C and will remain stable in the unopened vial until the expiration date. After opening, the contents should be used according to the instrument manufacturer's instructions and immediately returned to 2-8°C.

When used to monitor the precision of laboratory testing procedures for BNP and NT-proBNP assays, Control FD BNP has a reconstituted stability of up to 5 days for NT-ProBNP and 2 days for BNP under the proper storage conditions. Leaving the vial uncapped, or prolonging its time at room temperature, will void this open vial stability claim. Make sure the contents of the vial are well mixed before use.

PROCEDURE

Follow the manufacturer's instructions provided for BNP and NT-proBNP procedures. Verify that the lot number on the vial matches the assay sheet. To avoid evaporation, do not leave the vial uncapped. Controls should be run:

1. daily, in conjunction with patient samples.
2. as recommended by the instrument manufacturer.
3. as required by the relevant regulatory agency.

INSTRUCTIONS FOR USE

1. Remove a vial from the package.
2. Using a pipette, reconstitute the product with 1 mL of deionized water.
3. Allow the vial to sit at room temperature for 5 minutes.
4. Occasionally swirl for 15 minutes, or until all visible material is dissolved. Do not shake. Do not mix mechanically. Avoid getting any undissolved material on the sides of the vial or the stopper.
5. When all visible solid material is dissolved, invert several times to dissolve any material on the stopper.
6. Swirl occasionally for at least 5 minutes.
7. Use immediately or return to 2-8°C.
8. The vial should remain stored at 2-8°C at all times. If additional sampling is necessary, the time outside of 2-8°C storage should be minimized.

EXPECTED VALUES

The performance range for each level, based on data by combining estimates of assay variance as determined by participating laboratories using approved FDA instruments and reagents, is provided below. Average values obtained in the laboratory should fall within the performance range although the recovery may not be identical with the mean value listed. Variation between labs will be greater than the precision for any one instrument. Accuracy and precision depend on differences in equipment, reagents, supplies and techniques. Therefore, a lab must establish its own acceptable target values and ranges.

LIMITATIONS OF THE PROCEDURE

Variations in instruments and in the temperature of the testing material may result in accuracy differences. If the liquid in the vial becomes frozen, discard and use another vial, as the results will not be valid.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

ORDERING INFORMATION

PRODUCT NUMBER	PRODUCT DESCRIPTION	PRODUCT PACKAGING
K023M-6	Control FD BNP	6 x 1 mL

Distributed by Audit® MicroControls™, Inc. - U.S. customers only please call (866) 252-8348 or www.auditmicro.com

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