
Control LQ HbA1c

REF K067M-8

8 x 1.0 mL

LOT 07194

 2027-12-02

Aalto Scientific Ltd
230 Technology Pkwy
Eatonton, GA 31024
USA

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ENGLISH

INTENDED USE

Control LQ HbA1c is a reference control consisting of human blood based solutions. It is intended to simulate human patient whole blood samples for the purpose of monitoring the precision of laboratory testing procedures for the following assays: Glycohemoglobin A1c. When used for quality control purposes, it is recommended that each laboratory establish its own means and acceptable ranges and use the values provided only as guides.

SUMMARY AND PRINCIPLE

Good laboratory practices require that stable reference materials be used to verify the accuracy and precision of testing methods and techniques. Control LQ HbA1c may be used as one would use human whole blood to obtain the stated values.

REAGENTS

Control LQ HbA1c is an in vitro diagnostic control composed of a buffered bacteriostatic and fungistatic human blood matrix. Control LQ HbA1c is a liquid product manufactured according to standard quality control procedures. The manufacturer guarantees stability and consistency of this product.

WARNINGS AND PRECAUTIONS

Because this product is of human origin, it should be handled as though capable of transmitting infectious diseases. Each serum, plasma or whole blood donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled at the Centers for Disease Control's Biosafety Level 2.

Control LQ HbA1c is intended solely for in vitro diagnostic use for the purpose described on the labeling. Audit® MicroControls™, Inc. shall not be liable for any unclaimed damages arising from any other usage.

STORAGE AND STABILITY

Control LQ HbA1c is stored at 2-8°C and will remain stable in the unopened tube until the expiration date. After opening, the contents should be used according to the instrument manufacturer's instructions and immediately returned to 2-8°C.

When used to monitor the precision of laboratory testing procedures for Glycohemoglobin A1c assays, Control LQ HbA1c has an open vial stability of up to 50 days under the proper storage conditions. Leaving the tube uncapped, or prolonging its time at room temperature, will void this open vial stability claim. Make sure the contents of the tube are well mixed before use.

MATERIALS PROVIDED

Control LQ HbA1c, 8 x 1.0 mL

PROCEDURE

Follow the manufacturer's instructions provided for Glycohemoglobin A1c procedures. Verify that the lot number on the tube matches the assay sheet. To avoid evaporation, do not leave the tube uncapped. Controls should be run:

1. daily, in conjunction with patient samples.
2. as recommended by the instrument manufacturer.
3. as required by the relevant regulatory agency.

INSTRUCTIONS FOR USE

1. Remove the tubes from the package and mix by gentle inversion. Do not shake. Do not mix mechanically.
2. Gently swirl to assure complete mixing of the contents.
3. Unscrew the caps and remove them from the tubes. (Caps may be left on if analyzer is capable of puncturing a vacutainer top.)
4. Run tests according to analyzer manufacturer instructions.
5. After testing is completed, replace the caps, then store at 2-8°C for 50 days.

EXPECTED VALUES

The performance range for each level, based on data by combining estimates of assay variance as determined by participating laboratories using approved FDA instruments and reagents, is provided below. Average values obtained in the laboratory should fall within the performance range although the recovery may not be identical with the mean value listed. Variation between labs will be greater than the precision for any one instrument. Accuracy and precision depend on differences in equipment, reagents, supplies and techniques. Therefore, a lab must establish its own acceptable target values and ranges.

U.S. customers only - Once each level of control is tested, raw data may be entered via the AUDITOR™ QC Program at www.auditmicro.com. An on-line graph showing actual values Levey-Jennings graphs and charts for each analyte is then available to print. Call (866) 25-AUDIT for more information.

LIMITATIONS OF THE PROCEDURE

Variations in instruments and in the temperature of the testing material may result in accuracy differences.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

ORDERING INFORMATION

PRODUCT NUMBER	PRODUCT DESCRIPTION	PRODUCT PACKAGING
K067M-8	Control LQ HbA1c	8 x 1.0 mL

Distributed by AUDIT® MicroControls™ Inc. - U.S. customers only please call (866) 252-8348 or www.auditmicro.com

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