
Control FD Procalcitonin

REF K071M-6

6 x 2 mL

LOT 07149

 2028-08-11

Aalto Scientific Ltd
230 Technology Parkway
Eatonton, GA 31024
USA

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ENGLISH

INTENDED USE

Control FD Procalcitonin is a reference control consisting of bovine serum albumin based solutions. It is intended to simulate human patient serum samples for the purpose of monitoring the precision of laboratory testing procedures for Procalcitonin. When used for quality control purposes, it is recommended that each laboratory establish its own means and acceptable ranges and use the values provided only as guides.

SUMMARY AND PRINCIPLE

Good laboratory practices require that stable reference materials be used to verify the accuracy and precision of testing methods and techniques. Control FD Procalcitonin may be used as one would use human serum to obtain the stated values.

REAGENTS

Control FD Procalcitonin is an in vitro diagnostic control composed of a buffered bacteriostatic and fungistatic bovine serum albumin matrix. Control FD Procalcitonin is a lyophilized product manufactured according to standard quality control procedures. The manufacturer guarantees stability and consistency of this product.

WARNINGS AND PRECAUTIONS

This product contains less than 0.1% sodium azide that may react with lead and copper plumbing to form potentially explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Control FD Procalcitonin is intended solely for in vitro diagnostic use for the purpose described on the labeling. Audit® MicroControls™, Inc. shall not be liable for any unclaimed damages arising from any other usage.

STORAGE AND STABILITY

Control FD Procalcitonin is stored at 2-8°C and will remain stable in the unopened vial until the expiration date. After opening, the contents should be used according to the instrument manufacturer's instructions and immediately returned to 2-8°C.

When used to monitor the precision of laboratory testing procedures for Procalcitonin assays, Control FD Procalcitonin has a reconstituted stability of up to 7 days under the proper storage conditions. Leaving the vial uncapped, or prolonging its time at room temperature, will void this open vial stability claim. Make sure the contents of the vial are well mixed before use.

MATERIALS PROVIDED

Control FD Procalcitonin, 6 x 2 mL

PROCEDURE

Follow the manufacturer's instructions provided for Procalcitonin procedures. Verify that the lot number on the vial matches the assay sheet. To avoid evaporation, do not leave the vial uncapped. Controls should be run:

1. daily, in conjunction with patient samples.
2. as recommended by the instrument manufacturer.
3. as required by the relevant regulatory agency.

INSTRUCTIONS FOR USE

1. Remove a vial from the package.
2. Using a pipette, reconstitute the product with 2 mL of deionized water.
3. Allow the vial to sit at room temperature for 5 minutes.
4. Occasionally swirl for 15 minutes, or until all visible material is dissolved. Do not shake. Do not mix mechanically. Avoid getting any undissolved material on the sides of the vial or the stopper.
5. When all visible solid material is dissolved, invert several times to dissolve any material on the stopper.
6. Swirl occasionally for at least 5 minutes.
7. Use immediately or return to 2-8°C.
8. The vial should remain stored at 2-8°C at all times. If additional sampling is necessary, the time outside of 2-8°C storage should be minimized.

EXPECTED VALUES

The performance range for each level, based on data by combining estimates of assay variance as determined by participating laboratories using approved FDA instruments and reagents, is provided below. Average values obtained in the laboratory should fall within the performance range although the recovery may not be identical with the mean value listed. Variation between labs will be greater than the precision for any one instrument. Accuracy and precision depend on differences in equipment, reagents, supplies and techniques. Therefore, a lab must establish its own acceptable target values and ranges.

LIMITATIONS OF THE PROCEDURE

Variations in instruments and in the temperature of the testing material may result in accuracy differences. Make sure that the vial is brought to room temperature before testing. If the liquid in the vial becomes frozen, discard and use another vial, as the results will not be valid.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

ORDERING INFORMATION

PRODUCT NUMBER	PRODUCT DESCRIPTION	PRODUCT PACKAGING
K071M-6	Control FD Procalcitonin	6 x 2 mL


Distributed by Audit® MicroControls™, Inc. - U.S. customers only please call (866) 252-8348 or www.auditmicro.com

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