
Control LQ RF/CRP

REF K073M-6

6 x 2 mL

LOT 06930B

 2027-07-17

Aalto Scientific Ltd
230 Technology Park Way
Eatonton, GA 31024
USA

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ENGLISH

INTENDED USE

Control LQ RF/CRP is an assayed, stable, ready-to-use liquid, tri-level reference control for use with assays designed to quantitate Rheumatoid Factor (RF) and C-Reactive Protein (CRP). It is intended to simulate human patient samples for the purpose of monitoring the precision of laboratory testing procedures for RF and CRP assays. When used for quality control purposes, it is recommended that each laboratory establish its own means and acceptable ranges and use the values provided only as guides.

This product is intended for use with quantitative assays on the indicated analyzer provided in the labeling. The Control LQ RF/CRP should not be used for calibration or standardization of the RF and CRP assay. Control LQ RF/CRP is for In Vitro Diagnostic use only.

SUMMARY AND PRINCIPLE

Good laboratory practices require that stable reference materials be used to verify the accuracy and precision of testing methods and techniques. Control LQ RF/CRP may be used as one would use human serum to obtain the stated RF and CRP values. This control will assist in the evaluation of proper performance of RF and CRP assays.

REAGENTS

Control LQ RF/CRP is prepared using human serum with purified RF and CRP added to achieve the desired concentration levels, and 0.1% sodium azide as a preservative. Control LQ RF/CRP is a ready to use liquid product manufactured according to standard quality control procedures. The manufacturer guarantees stability and consistency of this product.

WARNINGS AND PRECAUTIONS

Because this product is of human origin, it should be handled as though capable of transmitting infectious diseases. Each serum, plasma or whole blood donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled at the Centers for Disease Control's Biosafety Level 2.

This product contains less than 0.1% sodium azide that may react with lead and copper plumbing to form potentially explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Control LQ RF/CRP is intended solely for in vitro diagnostic use for the purpose described on the labeling. Audit® MicroControls™, Inc. shall not be liable for any unclaimed damages arising from any other usage.

MATERIALS PROVIDED

Control LQ RF/CRP, 6 x 2 mL

STORAGE AND STABILITY

Control LQ RF/CRP is stored at 2-8°C and will remain stable in the unopened vial until the expiration date. After opening, the contents should be used according to the instrument manufacturer's instructions and immediately returned to 2-8°C.

When used to monitor the precision of laboratory testing procedures for RF and CRP assays, Control LQ RF/CRP has an open vial stability of up to 30 days under the proper storage conditions. Leaving the vial uncapped, or prolonging its time at room temperature, will void this open vial stability claim. Make sure the contents of the vial are well mixed before use.

PROCEDURE

Follow the manufacturer's instructions provided for RF and CRP procedures. Verify that the lot number on the vial matches the assay sheet. To avoid evaporation, do not leave the vial uncapped. Controls should be run:

1. daily, in conjunction with patient samples.
2. as recommended by the instrument manufacturer.
3. as required by the relevant regulatory agency.

INSTRUCTIONS FOR USE

1. Remove a vial from the package and mix by gentle inversion. Do not shake. Do not mix mechanically.
2. Use immediately or return to 2-8°C.
3. The vial should remain stored at 2-8°C at all times. If additional sampling is necessary, the time outside of 2-8°C storage should be minimized.

EXPECTED VALUES

The performance range for each level, based on data by combining estimates of assay variance as determined by participating laboratories using approved FDA instruments and reagents, is provided below. Average values obtained in the laboratory should fall within the performance range although the recovery may not be identical with the mean value listed. Variation between labs will be greater than the precision for any one instrument. Accuracy and precision depend on differences in equipment, reagents, supplies and techniques. Therefore, a lab must establish its own acceptable target values and ranges.

LIMITATIONS OF THE PROCEDURE

If the liquid in the vial becomes frozen, discard and use another vial, as the results will not be valid.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

ORDERING INFORMATION

PRODUCT NUMBER	PRODUCT DESCRIPTION	PRODUCT PACKAGING
K073M-6	Control LQ RF/CRP	6 x 2 mL

Distributed by , Audit® MicroControls™ Inc. - U.S. customers only please call (866) 252-8348 or www.auditmicro.com

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